

EdemaWear® Sizing Guide

ARM or LEG

Size	Arm/Calf/Thigh Circumference	Site
Small (#0600001) Navy Stripe	Up to 18" (45cm) Stockinet Length: 22"	* Wrist to Shoulder * Foot to Knee
Medium (#0960001) Yellow Stripe	Up to 30" (75cm) Stockinet Length: 34"	* Wrist to Shoulder * Foot to Groin ** may be cut to fit for below knee use
Large (#120L001) Red Stripe	Up to 46" (115cm) Stockinet Length: 24"	* Knee to groin for obese individuals
X-Large (#160XL01) Aqua Stripe	Up to 60" (150cm) Stockinet Length: 24"	* Knee to groin for morbidly obese individuals

EdemaWear® LITE for Tender Skin

Small (#L060017) Purple Stripe	Up to 24" (60cm) Stockinet Length: 22"	* Wrist to Shoulder * Foot to Knee
Medium (#L096034) Orange Stripe	Up to 36" (90cm) Stockinet Length: 34"	* Wrist to Shoulder * Foot to Groin ** may be cut to fit for below knee use

EdemaWear® SHAPED

Small (#SH06018) Grey Stripe	Up to 24" (60cm) Stockinet Length: 30"	* Foot to Mid-Thigh
------------------------------------	---	---------------------

Made in the U.S.A.

Revised 020216

Directions for Use:

- * Roll the EdemaWear stockinet up in your hands like you would a stocking before putting it on.
- * EdemaWear can be placed on a swollen or painful extremity by placing a plastic bag on the foot and pulling the EdemaWear stockinet over and above the plastic bag. Bag is then removed.
- * Remove daily and reapply.
- * To prevent EdemaWear stockinet from rolling down or up on the limb, turn over the top of EdemaWear making a 3" cuff, or turn EdemaWear inside out.
- * Cut to needed length. EdemaWear will not unravel.

Recommended Care:

- * Machine wash EdemaWear in cold water.
- * Do not use bleach.
- * Hang dry. Do not put in dryer.

Materials:

- * Medical Grade Lycra Spandex
- * Nylon
- * **NOT MADE** with Natural Rubber Latex

HCPCS:

- * A6457 when used on an open wound

U.S. Patent Numbers
8,034,013B2; 8,641,653B2; 9,259,373
Patent Pending in Canada
WIPO/PCT# WO 2014/144790A1



Omaha, Nebraska
Tel: (402) 660-9834 Fax: (402) 505-4690
Email: info@EdemaWear.com
www.EdemaWear.com