

EdemaWear® Sizing Guide

ARM or LEG

Size	Arm/Calf/Thigh Circumference	Site
Small (#0600001) Navy Stripe	Up to 18" (45cm) <small>Stockinet Length: 22"</small>	* Wrist to Shoulder * Foot to Knee
Medium (#0960001) Yellow Stripe	Up to 30" (75cm) <small>Stockinet Length: 34"</small>	* Wrist to Shoulder * Foot to Groin ** may be cut to fit for below knee use
Large (#120L001) Red Stripe	Up to 46" (115cm) <small>Stockinet Length: 24"</small>	* Knee to groin for obese individuals
X-Large (#160XL01) Aqua Stripe	Up to 60" (150cm) <small>Stockinet Length: 24"</small>	* Knee to groin for morbidly obese individuals

EdemaWear® LITE for Tender Skin

Small (#L060017) Purple Stripe	Up to 24" (60cm) <small>Stockinet Length: 22"</small>	* Wrist to Shoulder * Foot to Knee
Medium (#L096034) Orange Stripe	Up to 36" (90cm) <small>Stockinet Length: 34"</small>	* Wrist to Shoulder * Foot to Groin ** may be cut to fit for below knee use

EdemaWear® SHAPED

Small (#SH06018) Grey Stripe	Up to 24" (60cm) <small>Stockinet Length: 30"</small>	* Foot to Mid-Thigh
------------------------------------	--	---------------------